



uxuA Antibody

Product Code	CSB-PA432685ZA01SUD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A9EVN4
Immunogen	Recombinant Sorangium cellulosum uxuA protein
Raised In	Rabbit
Species Reactivity	Sorangium cellulosum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mannonate dehydratase (EC 4.2.1.8) (D-mannonate hydro-lyase), uxuA
Immunogen Species	Sorangium cellulosum (strain So ce56) (Polyangium cellulosum (strain So ce56))
Lead Time	Made-to-order (12-14 weeks)
Target Names	uxuA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.