



pcrA Antibody

Product Code	CSB-PA349564ZA01GFM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56255
Immunogen	Recombinant Geobacillus stearothermophilus pcrA protein
Raised In	Rabbit
Species Reactivity	Geobacillus stearothermophilus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent DNA helicase PcrA (EC 3.6.4.12), pcrA
Immunogen Species	Geobacillus stearothermophilus (Bacillus stearothermophilus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	pcrA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.