



Levanase Antibody

Product Code	CSB-PA523248ZA01BRG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O31411
Immunogen	Recombinant Bacillus sp. Levanase protein
Raised In	Rabbit
Species Reactivity	Bacillus sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Levanase (EC 3.2.1.65) (2,6-beta-D-fructan fructanohydrolase) (Endo-levanase) (Fragment)
Immunogen Species	Bacillus sp. (strain L7)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.