



Mite group 2 allergen Lep d 2 Antibody

Product Code	CSB-PA304483ZA01LEA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80384
Immunogen	Recombinant Lepidoglyphus destructor Mite group 2 allergen Lep d 2 protein
Raised In	Rabbit
Species Reactivity	Lepidoglyphus destructor
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mite group 2 allergen Lep d 2 (Allergen Lep d 1) (Allergen Lep d I) (allergen Lep d 2)
Immunogen Species	Lepidoglyphus destructor (Storage mite) (Glycyphagus destructor)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.