



BLOT5 Antibody

Product Code	CSB-PA527280ZA01BTN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O96870
Immunogen	Recombinant Blomia tropicalis BLOT5 protein
Raised In	Rabbit
Species Reactivity	Blomia tropicalis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mite allergen Blo t 5 (allergen Blo t 5), BLOT5
Immunogen Species	Blomia tropicalis (Mite)
Lead Time	Made-to-order (12-14 weeks)
Target Names	BLOT5
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.