



DERF2 Antibody

Product Code	CSB-PA312422ZA01DCO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00855
Immunogen	Recombinant Dermatophagoides farinae DERF2 protein
Raised In	Rabbit
Species Reactivity	Dermatophagoides farinae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mite group 2 allergen Der f 2 (Allergen Der f II) (allergen Der f 2), DERF2
Immunogen Species	Dermatophagoides farinae (American house dust mite)
Lead Time	Made-to-order (12-14 weeks)
Target Names	DERF2
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.