




DERP2 Antibody

Product Code	CSB-PA343403ZA01DJH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49278
Immunogen	Recombinant Dermatophagoides pteronyssinus Mite group 2 allergen Der p 2 (DERP2) (18-146AA)
Raised In	Rabbit
Species Reactivity	Dermatophagoides pteronyssinus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Mite group 2 allergen Der p 2 (Allergen Der p II) (DPX) (allergen Der p 2), DERP2
Immunogen Species	Dermatophagoides pteronyssinus (European house dust mite)
Target Names	DERP2

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: DERP2 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 30.1 kDa

Observed band size: 38 kDa

Product Modify Unmodified

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.