



IptD Antibody

Product Code	CSB-PA845655ZA01SWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8Z9J6
Immunogen	Recombinant Salmonella typhi IptD protein
Raised In	Rabbit
Species Reactivity	Salmonella typhi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	LPS-assembly protein LptD, IptD, imp ostA
Immunogen Species	Salmonella typhi
Lead Time	Made-to-order (12-14 weeks)
Target Names	IptD
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.