




DERP1 Antibody

Product Code	CSB-PA357315ZA01DJH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08176
Immunogen	Recombinant Dermatophagoides pteronyssinus Peptidase 1 (DERP1) (99-320AA)
Raised In	Rabbit
Species Reactivity	Dermatophagoides pteronyssinus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidase 1 (EC 3.4.22.65) (Allergen Der p I) (Major mite fecal allergen Der p 1) (allergen Der p 1), DERP1
Immunogen Species	Dermatophagoides pteronyssinus (European house dust mite)
Target Names	DERP1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: DERP1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41 kDa

Observed band size: 41 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.