



# Dsg3 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA007205ZA01MO   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| <b>Uniprot No.</b>         | O35902   |
| <b>Immunogen</b>           | Recombinant Mus musculus Dsg3 protein  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Mus musculus   |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)                                     |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Desmoglein-3 (130 kDa pemphigus vulgaris antigen homolog), Dsg3                  |
| <b>Immunogen Species</b>   | Mus musculus (Mouse)   |
| <b>Lead Time</b>           | Made-to-order (12-14 weeks)  |
| <b>Target Names</b>        | Dsg3   |
| <b>Product Modify</b>      | Unmodified   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.      |