



OLIG2 Antibody

Product Code	CSB-PA623814ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13516
Immunogen	Recombinant Homo sapiens OLIG2 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oligodendrocyte transcription factor 2 (Oligo2) (Class B basic helix-loop-helix protein 1) (bHLHb1) (Class E basic helix-loop-helix protein 19) (bHLHe19) (Protein kinase C-binding protein 2) (Protein kinase C-binding protein RACK17), OLIG2, BHLHB1 BHLHE19 PRKCBP2 RACK17
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	OLIG2
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.