



At3g57050 Antibody

Product Code	CSB-PA345319ZA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53780
Immunogen	Recombinant Arabidopsis thaliana At3g57050 protein
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cystathionine beta-lyase, chloroplastic (CBL) (EC 4.4.1.8) (Beta-cystathionase) (Cysteine lyase)
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	At3g57050
Product Modify	Unmodified