



EURM2 Antibody

Product Code	CSB-PA874707ZA01EPP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9TZZ2
Immunogen	Recombinant Euroglyphus maynei Mite group 2 allergen Eur m 2 (EURM2) (17-145AA)
Raised In	Rabbit
Species Reactivity	Euroglyphus maynei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mite group 2 allergen Eur m 2 (allergen Eur m 2), EURM2
Immunogen Species	Euroglyphus maynei (Mayne's house dust mite)
Lead Time	Made-to-order
Target Names	EURM2
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.