



Toxin MIT1 Antibody

Product Code	CSB-PA326235ZA01DBI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25687
Immunogen	Recombinant Dendroaspis polylepis polylepis Toxin MIT1 protein
Raised In	Rabbit
Species Reactivity	Dendroaspis polylepis polylepis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toxin MIT1 (MIT 1) (Black mamba intestinal toxin 1) (Black mamba venom protein A)
Immunogen Species	Dendroaspis polylepis polylepis (Black mamba)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.