



Glutathione S-transferase Antibody

Product Code	CSB-PA522732ZA01BTL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O18598
Immunogen	Recombinant Blattella germanica Glutathione S-transferase protein
Raised In	Rabbit
Species Reactivity	Blattella germanica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase (EC 2.5.1.18) (GST class-sigma) (Major allergen Bla g 5) (allergen Bla g 5)
Immunogen Species	Blattella germanica (German cockroach) (Blatta germanica)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.