



DDP1 Antibody

Product Code	CSB-PA860334ZA01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99321
Immunogen	Recombinant Saccharomyces cerevisiae DDP1 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diphosphoinositol polyphosphate phosphohydrolase DDP1 (EC 3.6.1.52) (Diadenosine 5',5'''-P1,P6-hexaphosphate hydrolase) (Ap6A hydrolase) (Diadenosine and diphosphoinositol polyphosphate phosphohydrolase 1) (Diadenosine hexaphosphate hydrolase (AMP-forming)) (EC 3.6.1.60), DDP1
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Lead Time	Made-to-order (12-14 weeks)
Target Names	DDP1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.