



MTHFD2 Antibody

Product Code	CSB-PA606691ZA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0P5C2
Immunogen	Recombinant Bos taurus MTHFD2 protein
Raised In	Rabbit
Species Reactivity	Bos taurus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bifunctional methylenetetrahydrofolate dehydrogenase/cyclohydrolase, mitochondrial [Includes: NAD-dependent methylenetetrahydrofolate dehydrogenase (EC 1.5.1.15); Methenyltetrahydrofolate cyclohydrolase (EC 3.5.4.9)], MTHFD2
Immunogen Species	Bos taurus (Bovine)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	MTHFD2
Product Modify	Unmodified