



MTHFD1 Antibody

Product Code	CSB-PA732682ZA01PYX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5R8P0
Immunogen	Recombinant Pongo abelii MTHFD1 protein
Raised In	Rabbit
Species Reactivity	Pongo abelii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-1-tetrahydrofolate synthase, cytoplasmic (C1-THF synthase) [Includes: Methylenetetrahydrofolate dehydrogenase (EC 1.5.1.5); Methenyltetrahydrofolate cyclohydrolase (EC 3.5.4.9); Formyltetrahydrofolate synthetase (EC 6.3.4.3)], MTHFD1
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Lead Time	Made-to-order (12-14 weeks)
Target Names	MTHFD1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.