



tmk Antibody

Product Code	CSB-PA007226ZA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WKE0
Immunogen	Recombinant Mycobacterium tuberculosis tmk protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thymidylate kinase (EC 2.7.4.9) (Thymidine monophosphate kinase) (dTMP kinase) (TMPK), tmk
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Lead Time	Made-to-order (12-14 weeks)
Target Names	tmk
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.