



Venom serine carboxypeptidase Antibody

Product Code	CSB-PA509855ZA01DNK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C9WMM5
Immunogen	Recombinant Apis mellifera Venom serine carboxypeptidase protein
Raised In	Rabbit
Species Reactivity	Apis mellifera
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Venom serine carboxypeptidase (EC 3.4.16.5) (allergen Api m 9)
Immunogen Species	Apis mellifera (Honeybee)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.