



Trx-2 Antibody

Product Code	CSB-PA892681ZA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9V429
Immunogen	Recombinant Drosophila melanogaster Trx-2 protein
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin-2 (DmTrx-2), Trx-2
Immunogen Species	Drosophila melanogaster (Fruit fly)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Trx-2
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.