



SH-1 Antibody

Product Code	CSB-PA361346ZA01ZAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04712
Immunogen	Recombinant Zea mays SH-1 protein
Raised In	Rabbit
Species Reactivity	Zea mays
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sucrose synthase 1 (EC 2.4.1.13) (Shrunken-1) (Sucrose-UDP glucosyltransferase 1), SH-1
Immunogen Species	Zea mays (Maize)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	SH-1
Product Modify	Unmodified