



EREG Antibody

Product Code	CSB-PA007779ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14944
Immunogen	Recombinant Homo sapiens EREG protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proepiregulin [Cleaved into: Epiregulin (EPR)], EREG
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	EREG
Product Modify	Unmodified