



AFUA_2G03830 Antibody

Product Code	CSB-PA525299ZA01NGS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60024
Immunogen	Recombinant Neosartorya fumigata AFUA_2G03830 protein
Raised In	Rabbit
Species Reactivity	Neosartorya fumigata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Allergen Asp f 4 (allergen Asp f 4)
Immunogen Species	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	AFUA_2G03830
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.