



K3L Antibody

Product Code	CSB-PA321831ZA01VAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20639
Immunogen	Recombinant Vaccinia virus K3L protein
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein K3
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	K3L
Product Modify	Unmodified