



HSP70 Antibody

Product Code	CSB-PA302100ZA01AZV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78983
Immunogen	Recombinant Alternaria alternata HSP70 protein
Raised In	Rabbit
Species Reactivity	Alternaria alternata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock 70 kDa protein (allergen Alt a 3) (Fragment), HSP70, ALTA3
Immunogen Species	Alternaria alternata (Alternaria rot fungus) (Torula alternata)
Lead Time	Made-to-order (12-14 weeks)
Target Names	HSP70
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.