



HA Antibody

Product Code	CSB-PA721113ZA01IFG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q67282
Immunogen	Recombinant Influenza A virus HA protein
Raised In	Rabbit
Species Reactivity	Influenza A virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemagglutinin [Cleaved into: Hemagglutinin HA1 chain; Hemagglutinin HA2 chain], HA
Immunogen Species	Influenza A virus (strain A/Mallard/New York/6750/1978 H2N2)
Lead Time	Made-to-order (12-14 weeks)
Target Names	HA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.