



HSP18.1 Antibody

Product Code	CSB-PA322513ZA01EWE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19243
Immunogen	Recombinant Pisum sativum HSP18.1 protein
Raised In	Rabbit
Species Reactivity	Pisum sativum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	18.1 kDa class I heat shock protein (HSP 18.1), HSP18.1
Immunogen Species	Pisum sativum (Garden pea)
Lead Time	Made-to-order (12-14 weeks)
Target Names	HSP18.1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.