



clpC Antibody

Product Code	CSB-PA857566ZA01SKX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99W78
Immunogen	Recombinant Staphylococcus aureus clpC protein
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent Clp protease ATP-binding subunit ClpC, clpC
Immunogen Species	Staphylococcus aureus (strain Mu50 / ATCC 700699)
Lead Time	Made-to-order (12-14 weeks)
Target Names	clpC
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.