



Defb1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA006662ZA01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P56386 |
| Immunogen | Recombinant Mus musculus Defb1 protein |
| Raised In | Rabbit |
| Species Reactivity | Mus musculus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Beta-defensin 1 (BD-1) (mBD-1) (Defensin, beta 1), Defb1 |
| Immunogen Species | Mus musculus (Mouse) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | Defb1 |
| Product Modify | Unmodified |