



IFNG Antibody

Product Code	CSB-PA011050ZA01MQG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35735
Immunogen	Recombinant Marmota monax IFNG protein
Raised In	Rabbit
Species Reactivity	Marmota monax
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon gamma (IFN-gamma), IFNG
Immunogen Species	Marmota monax (Woodchuck)
Lead Time	Made-to-order (12-14 weeks)
Target Names	IFNG
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.