



Tachylectin-2 Antibody

Product Code	CSB-PA634622ZA01TAH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q27084
Immunogen	Recombinant Tachypleus tridentatus Tachylectin-2 protein
Raised In	Rabbit
Species Reactivity	Tachypleus tridentatus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tachylectin-2 (Lectin L10c)
Immunogen Species	Tachypleus tridentatus (Japanese horseshoe crab)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified