



GPX4 Antibody

Product Code	CSB-PA670014ZA01EXP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4AEH2
Immunogen	Recombinant ongo pygmaeus GPX4 protein
Raised In	Rabbit
Species Reactivity	ongo pygmaeus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipid hydroperoxide glutathione peroxidase (PHGPx) (EC 1.11.1.12) (Glutathione peroxidase 4) (GPx-4) (GSHPx-4), GPX4
Immunogen Species	ongo pygmaeus (Bornean orangutan)
Lead Time	Made-to-order (12-14 weeks)
Target Names	GPX4
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.