



VACWR034 Antibody

Product Code	CSB-PA323585ZA01VAI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18378
Immunogen	Recombinant Vaccinia virus VACWR034 protein
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein K3 (Protein K2)
Immunogen Species	Vaccinia virus (strain Western Reserve) (VACV) (Vaccinia virus (strain WR))
Lead Time	Made-to-order (12-14 weeks)
Target Names	VACWR034
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.