



Tropomyosin Antibody

Product Code	CSB-PA631529ZA01DCO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q23939
Immunogen	Recombinant Dermatophagoides farinae Tropomyosin protein
Raised In	Rabbit
Species Reactivity	Dermatophagoides farinae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tropomyosin (Mag44) (allergen Der f 10)
Immunogen Species	Dermatophagoides farinae (American house dust mite)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.