



PHLPI Antibody

Product Code	CSB-PA331772ZA01EUQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43213
Immunogen	Recombinant Phleum pratense PHLPI protein
Raised In	Rabbit
Species Reactivity	Phleum pratense
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pollen allergen Phl p 1 (Allergen Phl p I) (allergen Phl p 1), PHLPI
Immunogen Species	Phleum pratense (Common timothy)
Lead Time	Made-to-order (12-14 weeks)
Target Names	PHLPI
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.