



Cyp17a1 Antibody

Product Code	CSB-PA006392ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11715
Immunogen	Recombinant Rattus norvegicus Cyp17a1 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Steroid 17-alpha-hydroxylase/17,20 lyase (EC 1.14.14.19) (EC 1.14.14.32) (17-alpha-hydroxyprogesterone aldolase) (CYPXVII) (Cytochrome P450 17A1) (Cytochrome P450-C17) (Cytochrome P450c17), Cyp17a1, Cyp17
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Cyp17a1
Product Modify	Unmodified