



LYZ Antibody

Product Code	CSB-PA013283ZA01LIQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24364
Immunogen	Recombinant Lophura leucomelanos LYZ protein
Raised In	Rabbit
Species Reactivity	Lophura leucomelanos
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysozyme C (EC 3.2.1.17) (1,4-beta-N-acetylmuramidase), LYZ
Immunogen Species	Lophura leucomelanos (Kalij pheasant)
Lead Time	Made-to-order (12-14 weeks)
Target Names	LYZ
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.