



# Tdo2 Antibody

<b>Product Code</b>	CSB-PA023351ZA01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P21643
<b>Immunogen</b>	Recombinant Rattus norvegicus Tdo2 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rattus norvegicus
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Tryptophan 2,3-dioxygenase (TDO) (EC 1.13.11.11) (Tryptamin 2,3-dioxygenase) (Tryptophan oxygenase) (TO) (TRPO) (Tryptophan pyrrolase) (Tryptophanase), Tdo2, Tdo
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	Tdo2
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.