



2S albumin Antibody

Product Code	CSB-PA20169ZA01CBAD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D0PWG2
Immunogen	Recombinant Corylus avellana 2S albumin protein
Raised In	Rabbit
Species Reactivity	Corylus avellana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	2S albumin
Immunogen Species	Corylus avellana (European hazel) (Corylus maxima)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.