



ypr-10 Antibody

Product Code	CSB-PA20179ZA01AXF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D1YSM5
Immunogen	Recombinant Actinidia deliciosa ypr-10 protein
Raised In	Rabbit
Species Reactivity	Actinidia deliciosa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bet v 1 related allergen (ypr-10)
Immunogen Species	Actinidia deliciosa (Kiwi)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ypr-10
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.