



GP Antibody

Product Code	CSB-PA311249ZA01ZAT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P87666
Immunogen	Recombinant Zaire ebolavirus GP protein
Raised In	Rabbit
Species Reactivity	Zaire ebolavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Envelope glycoprotein (GP1,2) (GP) [Cleaved into: GP1; GP2; GP2-delta], GP
Immunogen Species	Zaire ebolavirus (strain Kikwit-95) (ZEBOV) (Zaire Ebola virus)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	GP
Product Modify	Unmodified