



prsA Antibody

Product Code	CSB-PA702193ZA01FLB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5HET4
Immunogen	Recombinant Staphylococcus aureus prsA protein
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Foldase protein PrsA (EC 5.2.1.8), prsA
Immunogen Species	Staphylococcus aureus (strain COL)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	prsA
Product Modify	Unmodified