



N Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA741755ZA01NCT |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q77K03 |
| Immunogen | Recombinant Newcastle disease virus N protein |
| Raised In | Rabbit |
| Species Reactivity | Newcastle disease virus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Nucleoprotein (Nucleocapsid protein) (NP) (Protein N), N, NP |
| Immunogen Species | Newcastle disease virus (strain Chicken/United States/B1/48) (NDV) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | N |
| Product Modify | Unmodified |