



bopE Antibody

Product Code	CSB-PA714045ZA01BNS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q62B14
Immunogen	Recombinant Burkholderia mallei bopE protein
Raised In	Rabbit
Species Reactivity	Burkholderia mallei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide exchange factor BopE (Effector protein BopE), bopE
Immunogen Species	Burkholderia mallei (strain ATCC 23344)
Lead Time	Made-to-order (12-14 weeks)
Target Names	bopE
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.