



Calreticulin Antibody

Product Code	CSB-PA004458ZA01EKM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P83003
Immunogen	Recombinant Entamoeba histolytica Calreticulin protein
Raised In	Rabbit
Species Reactivity	Entamoeba histolytica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calreticulin (Fragment)
Immunogen Species	Entamoeba histolytica
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.