



rnfH Antibody

Product Code	CSB-PA481560ZA01DXM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B6IZH9
Immunogen	Recombinant Coxiella burnetii rnfH protein
Raised In	Rabbit
Species Reactivity	Coxiella burnetii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein RnfH, rnfH
Immunogen Species	Coxiella burnetii (strain CbuG_Q212) (Coxiella burnetii (strain Q212))
Lead Time	Made-to-order (12-14 weeks)
Target Names	rnfH
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.