



RC0497 Antibody

Product Code	CSB-PA849886ZA01RMS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92IC3
Immunogen	Recombinant Rickettsia conorii RC0497 protein
Raised In	Rabbit
Species Reactivity	Rickettsia conorii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative N-acetylmuramoyl-L-alanine amidase RC0497 (EC 3.5.1.28)
Immunogen Species	Rickettsia conorii (strain ATCC VR-613 / Malish 7)
Lead Time	Made-to-order (12-14 weeks)
Target Names	RC0497
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.