



JTB Antibody

Product Code	CSB-PA011970ZA01BMV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O77049
Immunogen	Recombinant Brugia malayi JTB protein
Raised In	Rabbit
Species Reactivity	Brugia malayi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Protein JTB, JTB
Immunogen Species	Brugia malayi (Filarial nematode worm)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	JTB
Product Modify	Unmodified