



PPIAL4A Antibody

Product Code	CSB-PA896522ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y536
Immunogen	Recombinant Homo sapiens PPIAL4A protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase A-like 4A (PPIase A-like 4A) (EC 5.2.1.8) (Chromosome one-amplified sequence 2) (COAS-2) (Cyclophilin homolog overexpressed in liver cancer), PPIAL4A, COAS2 PPIAL4B
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	PPIAL4A
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.